

Anti-Kir3DL3 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8N743

Item Number: CM0015637

Introduction

Anti-KIR3DL3 rabbit polyclonal antibody validated for WB, IHC-P, and ELISA applications. Recognizes human, mouse, and rat KIR3DL3, an NK cell inhibitory receptor. Ideal for immune cell and oncology research.

[Learn More](#)

Parameter	Value
Product Name	KIR3DL3 Rabbit pAb
Remarks/Alias	KIR44; KIRC1; CD158Z; KIR3DL7; KIR2DL5B; KIR3DL3
Species	Human
Gene ID (Human)	115653
Gene ID	115653
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB,IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N743
Protein Weight	45kDa
Shipping	Ice bag