

Anti-Kir3Ds1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14943

Item Number: CM0025098

Introduction

KIR3DS1 Rabbit pAb detects the NK cell receptor CD158E2 across human, mouse, and rat. Validated for Western blot, immunohistochemistry, and ELISA. Supports innate immunity and MHC class I receptor studies.

[Learn More](#)

Parameter	Value
Product Name	KIR3DS1 Rabbit pAb
Remarks/Alias	KIR-G1, NKAT10, CD158E2, NKAT-10, KIR-123FM, KIR3DS1
Species	Human
Gene ID (Human)	3813
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-340 of human KIR3DS1 (NP_001077008.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14943
Protein Weight	42kDa
Shipping	Ice bag