

Anti-Kirrel Polyclonal Antibody For Wb, Elisa - Q96J84

Item Number: CM0013109

Introduction

Rabbit polyclonal antibody against human KIRREL (NEPH1) validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for podocyte biology and glomerular filtration barrier research.

[Learn More](#)

Parameter	Value
Product Name	KIRREL Rabbit pAb
Remarks/Alias	NEPH1; KIRREL; NPHS23
Species	Human
GeneID (Human)	55243
GeneID	55243
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-200 of human KIRREL (NP_060710.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96J84
Protein Weight	84kDa
Shipping	Ice bag