

Anti-Kiss1 (Kisspeptin-1) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15726

Item Number: CM0028956



Introduction

Rabbit polyclonal antibody targeting human KISS1 (Kisspeptin-1, HH13). Validated for WB, IHC-P, and ELISA. Cross-reacts with human and mouse. Ideal for reproductive endocrinology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	KISS1 Rabbit pAb
Remarks/Alias	HH13; KISS-1; KISS1
Species	Human
Gene ID (Human)	3814
Gene ID	3814
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 40-138 of human KISS1 (NP_002247.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15726
Protein Weight	15kDa
Shipping	Ice bag