


Anti-Kiss1R Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q969F8

Item Number: CM0025252



Introduction

This KISS1R rabbit polyclonal antibody (43kDa) detects human, mouse, rat. Validated for WB, IHC-P, IF/ICC, ELISA. Ideal for reproductive biology and GPCR research.

[Learn More](#)

Parameter	Value
Product Name	KISS1R Rabbit pAb
Remarks/Alias	HH8; CPPB1; GPR54; AXOR12; KISS-1R; HOT7T175; KISS1R
Species	Human
GeneID	84634 (Human)
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human KISS1R (NP_115940.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q969F8
Protein Weight	43kDa
Shipping	Ice bag