


Anti-Kitlg Polyclonal Antibody For Wb, If/Icc, Elisa - P21583

Item Number: CM0030740



Introduction

Rabbit polyclonal antibody against human KITLG (SCF, MGF). Validated for WB, IF/ICC, ELISA; reacts with human, mouse, rat. For research in stem cell biology, hematopoiesis, and cancer signaling.

[Learn More](#)

Parameter	Value
Product Name	KITLG Rabbit pAb
Remarks/Alias	SF; MGF; SCF; SLF; DCUA; FPH2; FPHH; KL-1; Kitl; WS2F; SHEP7; DFNA69; KITLG
Species	Human
GeneID (Human)	4254
GeneID	4254
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21583
Protein Weight	31kDa
Shipping	Ice bag