


Anti-Kitlg Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P21583

Item Number: CM0019098



Introduction

Rabbit polyclonal antibody against human KITLG (SCF/c-Kit ligand). Validated in WB, IF/ICC, and ELISA with cross-reactivity to mouse. For hematopoietic, developmental, and oncology research.

[Learn More](#)

Parameter	Value
Product Name	KITLG Rabbit pAb
Remarks/Alias	SF, MGF, SCF, SLF, DCUA, FPH2, FPHH, KL-1, Kitl, WS2F, SHEP7, DFNA69, KITLG
Species	Human
Gene ID (Human)	4254
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-214 of human KITLG (NP_000890.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P21583
Protein Weight	31kDa
Shipping	Ice bag