

Anti-Klb Polyclonal Antibody For Wb, Elisa - Q86Z14

Item Number: CM0015857



Introduction

KLB Rabbit pAb, validated for WB and ELISA, detects human, mouse, and rat beta-klotho (KLB). Recombinant protein immunogen, 120kDa, suitable for studying FGF21 signaling and bile acid metabolism.

[Learn More](#)

Parameter	Value
Product Name	KLB Rabbit pAb
Remarks/Alias	BKL; KLB
Species	Human
Gene ID (Human)	152831
Gene ID	152831
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 720-990 of human KLB (NP_783864.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86Z14
Protein Weight	120kDa
Shipping	Ice bag