

Anti-Klc2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9H0B6

Item Number: CM0022455



Introduction

Rabbit polyclonal antibody targeting human KLC2 (kinesin light chain 2), cross-reactive with mouse and rat. Suitable for WB and ELISA, ideal for research on microtubule-based organelle transport and lysosome positioning in neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	KLC2 Rabbit pAb
Remarks/Alias	KLC2; SPOAN
Species	Human
Gene ID (Human)	64837
Gene ID	64837
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human KLC2 (NP_073733.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H0B6
Protein Weight	69kDa
Shipping	Ice bag