

Anti-Klf1 Polyclonal Antibody For Wb, Elisa - Q13351

Item Number: CM0014947



Introduction

KLF1 Rabbit Polyclonal Antibody (pAb) for Western Blot (WB) and ELISA applications. Cross-reacts with rat. Research use for erythrocyte development and transcription regulation.

[Learn More](#)

Parameter	Specification
Product Name	KLF1 Rabbit pAb
Remarks/Alias	EKLF; EKLF/KLF1; KLF1
Species	Human
GeneID (Human)	10661
GeneID	10661
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q13351
Protein Weight	38kDa
Shipping	Ice bag