

Anti-Klf1 Rabbit Polyclonal Antibody For Wb, Elisa - Q13351

Item Number: CM0022347



Introduction

KLF1 rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Reacts with human, mouse, and rat. Targets the erythroid transcription factor KLF1 (38 kDa), a key regulator of erythropoiesis and fetal-to-adult globin switching.

[Learn More](#)

Parameter	Value
Product Name	KLF1 Rabbit pAb
Remarks/Alias	KLF1; EKLF; CDAN4; HBFQTL6; INLU; Kruppel like factor 1
Species	Human
Gene ID	10661
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 211-311 of human KLF1 (NP_006554.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q13351
Protein Weight	38kDa
Shipping	Ice bag