

Anti-Klf10 Rabbit Monoclonal Antibody For Wb, Elisa - Q13118

Item Number: CM0006666



Introduction

Rabbit monoclonal antibody against human KLF10 (TIEG1), validated for Western blot and ELISA. Reacts with human and mouse samples. Suitable for circadian rhythm and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	KLF10 Rabbit mAb
Remarks/Alias	EGRA, TIEG, TIEG1, EGR-alpha, KLF10
Species	Human
Gene ID (Human)	7071
Gene ID	7071
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human KLF10 (Q13118)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13118
Protein Weight	53kDa
Shipping	Ice bag