

Anti-Klf11 Polyclonal Antibody For Wb, Elisa - O14901

Item Number: CM0023107

Introduction

Rabbit polyclonal antibody against human KLF11, validated for WB and ELISA. Reacts with mouse and rat. Suitable for transcription factor studies, TGF-beta signaling, and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	KLF11 Rabbit pAb
Remarks/Alias	FKLF; FKLF1; MODY7; TIEG2; Tieg3; KLF11
Species	Human
GeneID (human)	8462
GeneID	8462
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 13-107 of human KLF11 (NP_003588.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O14901
Protein Weight	55kDa
Shipping	Ice bag