

Anti-Klf12 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y4X4

Item Number: CM0023478

Introduction

KLF12 Rabbit pAb validated for WB, IHC-P, and ELISA; recognizes human, mouse, and rat KLF12 (AP2REP). Suitable for studying transcriptional repression and AP-2-alpha gene regulation.

[Learn More](#)

Parameter	Value
Product Name	KLF12 Rabbit pAb
Remarks/Alias	AP2REP; AP-2rep; HSPC122; KLF12
Species	Human
Gene ID	11278
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-310 of human KLF12 (NP_009180.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Y4X4
Protein Weight	44kDa
Shipping	Ice bag