

Anti-Klf13 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y2Y9

Item Number: CM0021738

Introduction

KLF13 Rabbit Polyclonal Antibody validated for WB and ELISA, targeting human, mouse, rat KLF13 transcription factor. Ideal for studying gene regulation, erythroid lineage, and T-cell activation.

[Learn More](#)

Parameter	Description
Product Name	KLF13 Rabbit pAb
Remarks/Alias	BTEB3; FKLF2; NSLP1; RFLAT1; RFLAT-1; KLF13
Species	Human
GeneID (Human)	51621
GeneID	51621
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human KLF13 (NP_057079.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2Y9
Protein Weight	31kDa
Shipping	Ice bag