

Anti-Klf14 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Td94

Item Number: CM0021391

Introduction

Rabbit polyclonal antibody targeting human KLF14 (BTEB5), validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Predicted molecular weight 33 kDa.

[Learn More](#)

Parameter	Value
Product Name	KLF14 Rabbit pAb
Remarks/Alias	BTEB5; KLF14
Species	Human
GeneID (Human)	136259
GeneID	136259
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 200-300 of human KLF14 (NP_619638.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TD94
Protein Weight	33kDa
Shipping	Ice bag