

## Anti-Klf15 Polyclonal Antibody For Wb, Elisa - Q9Uih9

Item Number: CM0030309



### Introduction

Rabbit polyclonal antibody against human KLF15, a transcriptional regulator involved in cardiac fibrosis and hypertrophy. Suitable for WB and ELISA. Cross-reacts with human. Ideal for cardiac and kidney research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | KLF15 Rabbit pAb                                                                                                 |
| Remarks/Alias    | KKLF; KLF15                                                                                                      |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 28999                                                                                                            |
| GeneID           | 28999                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human KLF15 (NP_054798.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human                                                                                                            |
| SWISS            | Q9UIH9                                                                                                           |
| Protein Weight   | 44kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |