

Anti-Klf15 Rabbit Polyclonal Antibody For If/Icc And Elisa - Q9Uih9

Item Number: CM0018078



Introduction

Rabbit polyclonal antibody against human KLF15 (Krueppel-like factor 15), validated for IF/ICC and ELISA. Suitable for studies on transcriptional regulation, cardiac fibrosis, and podocyte differentiation.

[Learn More](#)

Parameter	Value
Product Name	KLF15 Rabbit pAb
Remarks/Alias	KKLF; KLF15
Species	Human
GeneID (Human)	28999
GeneID	28999
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human KLF15 (NP_054798.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9UIH9
Protein Weight	44kDa
Shipping	Ice bag