

Anti-Klf4 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - O43474

Item Number: CM0002428



Introduction

KLF4 Rabbit Monoclonal Antibody validated for WB, IP, IF/ICC, IHC-P, ELISA. Detects human, mouse, rat KLF4 (O43474). Suitable for stem cell, differentiation, and epithelial barrier research.

[Learn More](#)

Parameter	Value
Product Name	KLF4 Rabbit mAb
Remarks/Alias	EZF; GKLf; KLF4
Species	Human
GeneID (human)	9314
GeneID	9314
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43474
Protein Weight	55 kDa
Shipping	Ice bag