

Anti-Klf4 Rabbit Polyclonal Antibody For Wb, Elisa - O43474

Item Number: CM0023907



Introduction

Rabbit polyclonal antibody targeting human KLF4, validated for Western blot and ELISA. Ideal for studying transcription regulation, stem cell maintenance, and epithelial differentiation.

[Learn More](#)

Parameter	Value
Product Name	KLF4 Rabbit pAb
Remarks/Alias	EZF; GKLF; KLF4
Species	Human
Gene ID (Human)	9314
Gene ID	9314
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-479 of human KLF4 (NP_004226.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O43474
Protein Weight	55kDa
Shipping	Ice bag