

Anti-Klf4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q60793

Item Number: CM0031005



Introduction

Rabbit polyclonal antibody targeting KLF4 transcription factor, validated for Western blot, immunohistochemistry-paraffin, and ELISA. Detects human, mouse, rat origins. Supports stem cell, epithelial biology, and developmental studies.

[Learn More](#)

Parameter	Value
Product Name	KLF4 Rabbit pAb
Remarks/Alias	EZF; Zie; Klf
Species	Mouse
GeneID (Human)	9314
GeneID	16600
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-115 of mouse KLF4 (NP_034767.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q60793
Protein Weight	52 kDa
Shipping	Ice bag