

Anti-Klf4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O43474

Item Number: CM0030968



Introduction

Rabbit anti-KLF4 polyclonal antibody validated for Western blot, immunohistochemistry, immunofluorescence, and ELISA. Reacts with human, mouse, rat. Suitable for developmental biology, stem cell, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	KLF4 Rabbit pAb
Remarks/Alias	EZF; GKLF; KLF4
Species	Human
Gene ID (Human)	9314
Gene ID	9314
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-115 of KLF4 (NP_001300981.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43474
Protein Weight	55kDa
Shipping	Ice bag