

Anti-Klf5 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13887

Item Number: CM0029857



Introduction

Krueppel-like factor 5 (KLF5) rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Detects human and mouse KLF5. Immunogen corresponds to the N-terminal region (aa 1-250). Suitable for transcription factor research. Uniprot accession Q13887.

[Learn More](#)

Parameter	Specification
Product Name	KLF5 Rabbit pAb
Remarks/Alias	CKLF, IKLF, BTEB2, KLF5
Species	Human
GeneID (Human)	688
Immunogen	Recombinant protein containing a sequence corresponding to amino acids 1-250 of human KLF5 (NP_001721.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13887
Protein Weight	51kDa
Shipping	Ice bag