

Anti-Klf7 Polyclonal Antibody For Wb, Elisa - 075840

Item Number: CM0027937

Introduction

CamelBio offers anti-KLF7 polyclonal antibody (Rabbit) validated for WB and ELISA. Detects human and mouse KLF7, a transcription factor involved in neuronal development, metabolism, and adipogenesis. Suitable for research use.

[Learn More](#)

Parameter	Value
Product Name	KLF7 Rabbit pAb
Remarks/Alias	UKLF; KLF7
Species	Human
Gene ID (Human)	8609
Gene ID	8609
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	075840
Protein Weight	33kDa
Shipping	Ice bag