

Anti-Klf9 Rabbit Polyclonal Antibody For Wb, Elisa - Q13886

Item Number: CM0028764

Introduction

Anti-KLF9 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, rat. Supports research on circadian rhythm, keratinocyte proliferation, and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	KLF9 Rabbit pAb
Remarks/Alias	BTEB; BTEB1; KLF9
Species	Human
GeneID (Human)	687
GeneID	687
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human KLF9 (NP_001197.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13886
Protein Weight	27kDa
Shipping	Ice bag