


Anti-Klhdc2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y2U9

Item Number: CM0016177



Introduction

KLHDC2 Rabbit Polyclonal Antibody validated for WB, IHC-P, and ELISA. Recognizes human and mouse KLHDC2, a substrate-recognition component of the CRL2 E3 ubiquitin ligase involved in protein quality control.

[Learn More](#)

Parameter	Specification
Product Name	KLHDC2 Rabbit pAb
Remarks/Alias	LCP; HCLP1; HCLP-1; KLHDC2
Species	Human
Gene ID (Human)	23588
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y2U9
Protein Weight	46kDa
Shipping	Ice bag