

Anti-Klhdc3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bq90

Item Number: CM0017227



Introduction

Rabbit polyclonal antibody against human KLHDC3, validated for Western blot, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Suitable for studying CRL2 ubiquitin ligase, C-terminal degron pathways, and meiosis.

[Learn More](#)

Parameter	Value
Product Name	KLHDC3 Rabbit pAb
Remarks/Alias	PEAS; dj20C7.3; KLHDC3
Species	Human
GeneID	116138
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-382 of human KLHDC3 (NP_476502.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BQ90
Protein Weight	43kDa
Shipping	Ice bag