

Anti-Klhl1 Polyclonal Antibody For Wb, Elisa - Q9Nr64

Item Number: CM0016107

Introduction

KLHL1 Rabbit pAb is a polyclonal antibody targeting human Kelch-like protein 1 (KLHL1), validated for Western blot and ELISA. KLHL1 is involved in actin cytoskeleton organization in brain cells.

[Learn More](#)

Parameter	Value
Product Name	KLHL1 Rabbit pAb
Remarks/Alias	MRP2; KLHL1
Species	Human
GeneID (Human)	57626
GeneID	57626
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human KLHL1 (NP_065917.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NR64
Protein Weight	83kDa
Shipping	Ice bag