


Anti-KLHL11 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nvr0

Item Number: CM0021481



Introduction

Rabbit polyclonal antibody targeting human KLHL11 (Kelch-like protein 11), a cullin-RING E3 ubiquitin ligase component; suitable for Western blotting and ELISA; cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	KLHL11 Rabbit pAb
Remarks/Alias	KLHL11
Species	Human
Gene ID (Human)	55175
Gene ID	55175
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 16-300 of human KLHL11 (NP_060613.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NVR0
Protein Weight	80kDa
Shipping	Ice bag