

Anti-KLHL15 Rabbit Polyclonal Antibody For Wb, Elisa - Q96M94

Item Number: CM0024077

Introduction

Anti-KLHL15 Rabbit Polyclonal Antibody validated for WB and ELISA, reactive with human, mouse, rat. Targets KLHL15, a substrate-specific adapter for CUL3 E3 ubiquitin ligase, involved in DNA damage repair and PP2A regulation.

[Learn More](#)

Parameter	Value
Product Name	KLHL15 Rabbit pAb
Remarks/Alias	XLID103; HEL-S-305; KLHL15
Species	Human
Gene ID	80311 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human KLHL15 (NP_085127.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96M94
Protein Weight	70kDa
Shipping	Ice bag