

Anti-KLHL3 Polyclonal Antibody For Wb, Elisa - Q9Uh77

Item Number: CM0013961



Introduction

KLHL3 Rabbit polyclonal antibody validated for Western blot and ELISA in human and mouse samples. Investigate renal ion transport, WNK kinase degradation, and ubiquitination pathways. Suitable for research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	KLHL3 Rabbit pAb
Remarks/Alias	PHA2D; KLHL3
Species	Human
Gene ID (Human)	26249
Gene ID	26249
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human KLHL3 (NP_059111.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UH77
Protein Weight	65kDa
Shipping	Ice bag