

Anti-KLHL42 Rabbit Polyclonal Antibody For IHC-P, Elisa - Q9P2K6

Item Number: CM0017071

Introduction

Anti-KLHL42 (Ctb9, KLHDC5) rabbit polyclonal antibody validated for IHC-P and ELISA applications. Reacts with human and mouse, targeting the Kelch-like protein 42 involved in E3 ubiquitin ligase complex and mitotic progression.

[Learn More](#)

Parameter	Value
Product Name	KLHL42 Rabbit pAb
Remarks/Alias	Ctb9; KLHDC5; KLHL42
Species	Human
GeneID (human)	57542
GeneID	57542
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 206-505 of human KLHL42 (NP_065833.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9P2K6
Protein Weight	57kDa
Shipping	Ice bag