

Anti-KLHL8 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9P2G9

Item Number: CM0016238

Introduction

KLHL8 Rabbit pAb: a polyclonal antibody raised in rabbit against human KLHL8, suitable for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat. Ideal for E3 ubiquitin ligase and RAPSN degradation studies.

[Learn More](#)

Parameter	Value
Product Name	KLHL8 Rabbit pAb
Remarks/Alias	KLHL8
Species	Human
GeneID (Human)	57563
GeneID	57563
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 271-620 of human KLHL8 (NP_065854.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P2G9
Protein Weight	69kDa
Shipping	Ice bag