

Anti-Klk1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P06870

Item Number: CM0022987



Introduction

Rabbit polyclonal antibody targeting human KLK1, validated for Western blot, IHC-P, and ELISA. Detects endogenous KLK1 in human, mouse, and rat samples. Ideal for kininogen processing studies.

[Learn More](#)

Parameter	Value
Product Name	KLK1 Rabbit pAb
Remarks/Alias	hK1; KLKR; Klk6; KLK1
Species	Human
GeneID (Human)	3816
GeneID	3816
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-262 of human KLK1 (NP_002248.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06870
Protein Weight	29kDa
Shipping	Ice bag