

Anti-Klk11 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Ubx7

Item Number: CM0026886

Introduction

Rabbit polyclonal antibody targeting human KLK11 (kallikrein-11, hippostasin). Validated for WB, IHC-P, ELISA; cross-reacts with human and mouse. Suitable for protease research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	KLK11 Rabbit pAb
Remarks/Alias	TLSP; PRSS20; KLK11
Species	Human
GeneID (Human)	11012
GeneID	11012
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UBX7
Protein Weight	31kDa
Shipping	Ice bag