

Anti-Klk15 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H2R5

Item Number: CM0024072

Introduction

KLK15 rabbit polyclonal antibody validated for WB and ELISA. Detects human KLK15 with cross-reactivity to mouse and rat. Ideal for protease research and biomarker studies.

[Learn More](#)

Parameter	Value
Product Name	KLK15 Rabbit pAb
Remarks/Alias	ACO; HSRNASPH; KLK15
Species	Human
GeneID	55554
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-180 of human KLK15 (NP_059979.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9H2R5
Protein Weight	28kDa
Shipping	Ice bag