


Anti-Klk3 Rabbit Polyclonal Antibody For Wb - P07288

Item Number: CM0022867



Introduction

Rabbit polyclonal anti-KLK3 (PSA) antibody, validated for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human and mouse prostate-specific antigen. Ideal for prostate cancer biomarker research.

[Learn More](#)

Parameter	Value
Product Name	KLK3 Rabbit pAb
Remarks/Alias	APS; PSA; hK3; KLK2A1; KLK3
Species	Human
GeneID (human)	354
GeneID	354
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 30-261 of human KLK3 (NP_001639.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P07288
Protein Weight	29kDa
Shipping	Ice bag