

Anti-Klk6 Polyclonal Antibody For Wb, Elisa - Q92876

Item Number: CM0016327



Introduction

KLK6 rabbit polyclonal antibody validated for Western blot (WB) and ELISA applications. Detects human kallikrein-6, a serine protease implicated in neurodegeneration and cancer metastasis.

[Learn More](#)

Parameter	Value
Product Name	KLK6 Rabbit pAb
Remarks/Alias	hK6; Bssp; Klk7; SP59; PRSS9; PRSS18; KLK6
Species	Human
Gene ID	5653
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human KLK6 (NP_002765.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q92876
Protein Weight	26kDa
Shipping	Ice bag