

Anti-Klk6 Rabbit Polyclonal Antibody For Wb, Elisa - Q92876

Item Number: CM0015958



Introduction

Rabbit polyclonal antibody against human Kallikrein-6 (KLK6, neurosin, PRSS9) validated for Western blot and ELISA. Detects KLK6 in human, mouse, and rat samples. Ideal for studies of neurodegeneration, cancer invasion, and protease biology.

[Learn More](#)

Parameter	Value
Product Name	KLK6 Rabbit pAb
Remarks/Alias	hK6; Bssp; Klk7; SP59; PRSS9; PRSS18; KLK6
Species	Human
GeneID (Human)	5653
GeneID	5653
Immunogen	Synthetic peptide, A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human KLK6 (NP_002765.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92876
Protein Weight	26kDa
Shipping	Ice bag