

## Anti-Klotho Polyclonal Antibody For Wb And Elisa - Q9Uef7

Item Number: CM0018685



### Introduction

KL Rabbit pAb detects Klotho in human, mouse, and rat by WB and ELISA. Suitable for studying anti-aging, phosphate homeostasis, and FGF23 signaling. Recombinant immunogen ensures specific detection.

[Learn More](#)

| Parameter        | Value                                                                                                          |
|------------------|----------------------------------------------------------------------------------------------------------------|
| Product Name     | KL Rabbit pAb                                                                                                  |
| Remarks/Alias    | KL; HFTC3; KL                                                                                                  |
| Species          | Human                                                                                                          |
| GeneID (Human)   | 9365                                                                                                           |
| GeneID           | 9365                                                                                                           |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human KL (NP_004786.2). |
| Source           | Rabbit                                                                                                         |
| Category         | Polyclonal Antibodies                                                                                          |
| Application      | WB, ELISA                                                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                                                              |
| SWISS            | Q9UEF7                                                                                                         |
| Protein Weight   | 116kDa                                                                                                         |
| Shipping         | Ice bag                                                                                                        |