


Anti-Klrc1 Polyclonal Antibody For Wb, Elisa - P26715

Item Number: CM0016832



Introduction

KLRC1 Rabbit pAb (P26715) validated for WB, ELISA. Recognizes human NKG2-A/NKG2-B, an immune inhibitory receptor on NK and T cells. Cross-reacts with mouse. For immuno-oncology and immune checkpoint research.

[Learn More](#)

Parameter	Value
Product Name	KLRC1 Rabbit pAb
Remarks/Alias	NKG2; NKG2A; CD159A; KLRC1
Species	Human
Gene ID (Human)	3821
Gene ID	3821
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 96-215 of human KLRC1 (NP_015567.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P26715
Protein Weight	26kDa
Shipping	Ice bag