

## Anti-Klrc4 Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - O43908

Item Number: CM0016794

### Introduction

High-quality rabbit polyclonal antibody targeting human KLRC4 (NKG2F), validated for western blot, IF/ICC, IF-P, and ELISA. Detects human, mouse, and rat NKG2-F protein. Ideal for NK cell receptor and MHC class I research.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | KLRC4 Rabbit pAb                                                                                                |
| Remarks/Alias    | NKG2F; NKG2-F; KLRC4                                                                                            |
| Species          | Human                                                                                                           |
| Gene ID (Human)  | 8302                                                                                                            |
| Gene ID          | 8302                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-55 of human KLRC4 (NP_038459.1) |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, IF/ICC, IF-P, ELISA                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                               |
| SWISS            | O43908                                                                                                          |
| Protein Weight   | 18kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |