

Anti-Klrd1 Polyclonal Antibody For Wb And Elisa - Q13241

Item Number: CM0022647



Introduction

High-quality rabbit polyclonal antibody targeting human KLRD1 (CD94), suitable for WB and ELISA with cross-reactivity to human and mouse. Ideal for investigating NK cell receptor biology and immune checkpoint mechanisms.

[Learn More](#)

Parameter	Value
Product Name	KLRD1 Rabbit pAb
Remarks/Alias	CD94; KLRD1
Species	Human
Gene ID (Human)	3824
Gene ID	3824
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-179 of human KLRD1 (NP_002253.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13241
Protein Weight	21kDa
Shipping	Ice bag