


Anti-Klrk1 Polyclonal Antibody For Wb, Elisa - P26718

Item Number: CM0028341



Introduction

Rabbit polyclonal antibody targeting human KLRK1/NKG2-D (SWISS: P26718), validated for WB and ELISA. Detects NKG2-D in human and mouse samples, ideal for immune surveillance and oncology research.

[Learn More](#)

Parameter	Value
Product Name	KLRK1 Rabbit pAb
Remarks/Alias	KLR; CD314; NKG2D; NKG2-D; D12S2489E; KLRK1
Species	Human
Gene ID (Human)	22914
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 77-216 of human KLRK1 (NP_031386.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P26718
Protein Weight	25kDa
Shipping	Ice bag