


Anti-Klrk1 Rabbit Polyclonal Antibody For If/Icc And Elisa - P26718

Item Number: CM0018331



Introduction

Rabbit polyclonal antibody targeting human KLRK1 (NKG2D), validated for IF/ICC and ELISA. Ideal for studying NK cell-mediated immunosurveillance and immune checkpoint research. Ships with ice pack.

[Learn More](#)

Parameter	Value
Product Name	KLRK1 Rabbit pAb
Remarks/Alias	KLR; CD314; NKG2D; NKG2-D; D12S2489E; KLRK1
Species	Human
Gene ID (Human)	22914
Gene ID	22914
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human KLRK1 (NP_031386.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P26718
Protein Weight	25kDa
Shipping	Ice bag