

Anti-Kmt1B / Suv39H2 Rabbit Monoclonal Antibody For Wb, Elisa - Q9H5I1

Item Number: CM0005208



Introduction

High-quality rabbit monoclonal antibody targeting human KMT1B/SUV39H2, a histone H3K9 methyltransferase critical for heterochromatin formation. Validated for WB and ELISA, suitable for epigenetic and cancer research.

[Learn More](#)

Parameter	Value
Product Name	KMT1B / SUV39H2 Rabbit mAb
Remarks/Alias	KMT1B; KMT1B / SUV39H2
Species	Human
Gene ID	79723 (Human)
Immunogen	Synthetic peptide corresponding to amino acids 311-410 of human KMT1B/SUV39H2 (UniProt: Q9H5I1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9H5I1
Protein Weight	47kDa
Shipping	Ice bag