


Anti-Kmt5A/Setd8 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nqr1

Item Number: CM0030528



Introduction

Rabbit polyclonal antibody targeting human KMT5A/SETD8, suitable for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for epigenetic and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	KMT5A/SETD8 Rabbit pAb
Remarks/Alias	SET8; SET07; SETD8; PR-Set7; PR/SET07; KMT5A/SETD8
Species	Human
GeneID (Human)	387893
GeneID	387893
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-393 of human KMT5A/SETD8 (Q9NQR1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQR1
Protein Weight	43kDa
Shipping	Ice bag