


Anti-Kmt5A/Setd8 Rabbit Monoclonal Antibody (Kd Validated) For Wb, Elisa - Q9Nqr1

Item Number: CM0008076



Introduction

KD-validated rabbit monoclonal antibody against KMT5A/SETD8 for Western blot and ELISA. Detects human and mouse KMT5A/SETD8, a histone H4K20 methyltransferase involved in gene silencing and cell cycle regulation. Ideal for epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] KMT5A/SETD8 Rabbit mAb
Remarks/Alias	SET8; SET07; SETD8; PR-Set7; PR/SET07; [KO/KD Validated] KMT5A/SETD8
Species	Human
GeneID	387893 (Human)
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human KMT5A/SETD8 (Q9NQR1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NQR1
Protein Weight	43kDa
Shipping	Ice bag