


Anti-Kmt5A/Setd8 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nqr1

Item Number: CM0029677



Introduction

KMT5A/SETD8 Rabbit Polyclonal Antibody validated for Western Blot, Immunofluorescence, and ELISA in human, mouse, rat. Targets the histone methyltransferase KMT5A/SETD8 (43kDa), a critical mitotic and transcriptional repressor via H4K20 monomethylation.

[Learn More](#)

Parameter	Value
Product Name	KMT5A/SETD8 Rabbit pAb
Remarks/Alias	SET8; SET07; SETD8; PR-Set7; PR/SET07; KMT5A/SETD8
Species	Human
Gene ID (Human)	387893
Gene ID	387893
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 73-352 of human KMT5A/SETD8 (NP_065115.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQR1
Protein Weight	43kDa
Shipping	Ice bag