

Anti-Kmt5C Rabbit Polyclonal Antibody For If-P And Elisa - Q86Y97

Item Number: CM0027692



Introduction

KMT5C Rabbit Polyclonal Antibody (SUV420H2) for IF-P and ELISA. Detects human histone methyltransferase KMT5C, a key enzyme for H4K20 di- and trimethylation, involved in transcriptional repression and DNA repair.

[Learn More](#)

Parameter	Value
Product Name	KMT5C Rabbit pAb
Remarks/Alias	SUV420H2; Suv4-20h2
Species	Human
GeneID (Human)	84787
GeneID	84787
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-248 of human KMT5C (NP_116090.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Human
SWISS	Q86Y97
Protein Weight	52kDa
Shipping	Ice bag